

For Research Use Only

SH3TC2 Polyclonal antibody

Catalog Number: 20067-1-AP



Basic Information

Catalog Number:

20067-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_024577

GeneID (NCBI):

79628

UNIPROT ID:

Q8TF17

Full Name:

SH3 domain and tetratricopeptide repeats 2

Calculated MW:

145 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

SH3TC2, also named as KIAA1985, contains 1 SH3 domain and TPR repeats. Defects in SH3TC2 are the cause of Charcot-Marie-Tooth disease type 4C (CMT4C) [MIM:601596]. The antibody can recognize all the isoforms except isoform4 of SH3TC2.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

