

For Research Use Only

SH3TC2 Polyclonal antibody

Catalog Number: 20067-1-AP



Basic Information

Catalog Number: 20067-1-AP	GenBank Accession Number: NM_024577	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 79628	
Source: Rabbit	Full Name: SH3 domain and tetrapeptide repeats 2	
Isotype: IgG	Calculated MW: 145 kDa	

Applications

Tested Applications:
IHC, ELISA

Species Specificity:
human, mouse, rat

Background Information

SH3TC2, also named as KIAA1985, contains 1 SH3 domain and TPR repeats. Defects in SH3TC2 are the cause of Charcot-Marie-Tooth disease type 4C (CMT4C) [MIM:601596]. The antibody can recognize all the isoforms except isoform 4 of SH3TC2.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data