

For Research Use Only

APOL5-Specific Polyclonal antibody

Catalog Number: 20068-1-AP



Basic Information

Catalog Number:

20068-1-AP

Size:

150ul , Concentration: 260 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_030642

GeneID (NCBI):

80831

UNIPROT ID:

Q9BWW9

Full Name:

apolipoprotein L, 5

Calculated MW:

47 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

APOL5 belongs to the apolipoprotein L family. The pathological expression of the apolipoproteins during pregnancy is implicated in fetal growth retardation, preeclampsia, and the onset of adult atherosclerosis [PubMed: 11374903]. APOL5 may affect the movement of lipids in the cytoplasm or allow the binding of lipids to organelles. This antibody is specific to APOL5.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

