

For Research Use Only

SCARF2-Specific Polyclonal antibody



Catalog Number: 20073-1-AP

Basic Information

Catalog Number: 20073-1-AP	GenBank Accession Number: NM_182895	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 213 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 91179	
Source: Rabbit	Full Name: scavenger receptor class F, member 2	
Isotype: IgG	Calculated MW: 92 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

SCARF2, also named as SREC2 and SREPCR, is probable adhesion protein, which mediates homophilic and heterophilic interactions. In contrast to SCARF1, it poorly mediates the binding and degradation of acetylated low density lipoprotein (Ac-LDL). This antibody is specific to SCARF2.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data