

For Research Use Only

# SCARF2-Specific Polyclonal antibody

Catalog Number: 20073-1-AP



## Basic Information

**Catalog Number:**

20073-1-AP

**Size:**

150ul , Concentration: 213 µg/ml by Bradford method using BSA as the standard;

**Source:**

Rabbit

**Isotype:**

IgG

**GenBank Accession Number:**

NM\_182895

**GeneID (NCBI):**

91179

**UNIPROT ID:**

Q96GP6

**Full Name:**

scavenger receptor class F, member 2

**Calculated MW:**

92 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

ELISA

**Species Specificity:**

human, mouse, rat

## Background Information

SCARF2, also named as SREC2 and SREPCR, is probable adhesion protein, which mediates homophilic and heterophilic interactions. In contrast to SCARF1, it poorly mediates the binding and degradation of acetylated low density lipoprotein (Ac-LDL). This antibody is specific to SCARF2.

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

