

For Research Use Only

ALPK2-Specific Polyclonal antibody

Catalog Number: 20079-1-AP



Basic Information

Catalog Number: 20079-1-AP	GenBank Accession Number: NM_052947	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 247 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 115701	
Source: Rabbit	Full Name: alpha-kinase 2	
Isotype: IgG	Calculated MW: 237 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

ALPK2, also named as HAK, belongs to the protein kinase superfamily, Alpha-type protein kinase family and ALPK subfamily. ALPK2 is a kinase that recognizes phosphorylation sites in which the surrounding peptides have an alpha-helical conformation. This antibody is specific to ALPK2.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data