

For Research Use Only

CSTL1-Specific Polyclonal antibody, PBS Only

Catalog Number: 20086-1-PBS



Basic Information

Catalog Number:

20086-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_138283

GeneID (NCBI):

128817

UNIPROT ID:

Q9H114

Full Name:

cystatin-like 1

Calculated MW:

17 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human, mouse

Background Information

CSTL1, also named as RCET11, belongs to the cystatin family. The antibody is specific to CSTL1.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data