

For Research Use Only

RIMS4-Specific Polyclonal antibody, PBS Only

Catalog Number:20091-1-PBS



Basic Information

Catalog Number:

20091-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_182970

GeneID (NCBI):

140730

UNIPROT ID:

Q9H426

Full Name:

regulating synaptic membrane
exocytosis 4

Calculated MW:

29 kDa

Observed MW:

29 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

RIMS4, also named as C20orf190, regulates synaptic membrane exocytosis. The antibody is specific to RIMS4.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20091-1-AP (RIMS4-Specific antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20091-1-PBS in a different storage buffer formulation.