

For Research Use Only

# LDLRAD2-Specific Polyclonal antibody



Catalog Number: 20095-1-AP

## Basic Information

|  |   |  |
|--|---|--|
| <b>Catalog Number:</b><br>20095-1-AP   | <b>GenBank Accession Number:</b><br>NM_001013693                                  | <b>Purification Method:</b><br>Antigen affinity purification |
| <b>Size:</b><br>150ul , Concentration: 280 µg/ml by Bradford method using BSA as the standard; | <b>GeneID (NCBI):</b><br>401944   |  |
| <b>Source:</b><br>Rabbit   | <b>Full Name:</b><br>low density lipoprotein receptor class A domain containing 2 |  |
| <b>Isotype:</b><br>IgG   | <b>Calculated MW:</b><br>29 kDa   |  |

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human

## Background Information

LDLRAD2 belongs to the LDLR family. This antibody is specific to LDLRAD2.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data