

For Research Use Only

NBPF4-Specific Polyclonal antibody, PBS Only

Catalog Number:20099-1-PBS



Basic Information

Catalog Number:

20099-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001143989

GeneID (NCBI):

148545

UNIPROT ID:

Q96M43

Full Name:

neuroblastoma breakpoint family,
member 4

Calculated MW:

72 kDa

Observed MW:

80-90 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

NBPF4 belongs to the NBPF family. This antibody is specific to NBPF4.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20099-1-AP (NBPF 4-Specific antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20099-1-PBS in a different storage buffer formulation.