

For Research Use Only

# ZSWIM1 Polyclonal antibody

Catalog Number: 20106-1-AP



---

## Basic Information

|   |  |  |
|---|--|--|
| <b>Catalog Number:</b><br>20106-1-AP        | <b>GenBank Accession Number:</b><br>BC001672             | <b>Purification Method:</b><br>Antigen affinity purification |
| <b>Size:</b><br>150ul                       | <b>GeneID (NCBI):</b><br>90204                           |  |
| <b>Source:</b><br>Rabbit                    | <b>Full Name:</b><br>zinc finger, SWIM-type containing 1 |  |
| <b>Isotype:</b><br>IgG                      | <b>Calculated MW:</b><br>485 aa, 55 kDa                  |  |
| <b>Immunogen Catalog Number:</b><br>AG13739 |  |  |

---

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human, mouse, rat

---

## Background Information

---

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

---

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data