

For Research Use Only

# ARHGEF40 Polyclonal antibody, PBS Only

Catalog Number: 20110-1-PBS



## Basic Information

**Catalog Number:**

20110-1-PBS

**Size:**

100ug, Concentration: 1 mg/ml by Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG13886

**GenBank Accession Number:**

BC000084

**GeneID (NCBI):**

55701

**UNIPROT ID:**

Q8TER5

**Full Name:**

hypothetical protein FLJ10357

**Calculated MW:**

414 aa, 46 kDa

**Observed MW:**

150-160 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

Indirect ELISA

**Species Specificity:**

human, mouse

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data