

For Research Use Only

# HNRNPUL2-Specific Polyclonal antibody, PBS Only

Catalog Number: 20143-1-PBS



## Basic Information

**Catalog Number:**

20143-1-PBS

**Size:**

100ug, Concentration: 1 mg/ml by Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**GenBank Accession Number:**

NM\_001079559

**GeneID (NCBI):**

221092

**UNIPROT ID:**

Q1KMD3

**Full Name:**

heterogeneous nuclear ribonucleoprotein U-like 2

**Calculated MW:**

85 kDa

**Observed MW:**

100-110 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

WB, IF/ICC, Indirect ELISA

**Species Specificity:**

human

## Background Information

HNRNPUL2 also named as SAF-A2. This antibody is specific to HNRNPUL2.

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20143-1-AP (HNRNPUL2-Specific antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20143-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using HNRNPUL2-Specific antibody (20143-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 20143-1-PBS in a different storage buffer formulation.