

For Research Use Only

FMNL1 Polyclonal antibody

Catalog Number: 20152-1-AP **1 Publications**



Basic Information

| | | |
|---|--|--|
| Catalog Number: 20152-1-AP | GenBank Accession Number: BC021906 | Purification Method: Antigen affinity purification |
| Size: 150ul | GeneID (NCBI): 752 | |
| Source: Rabbit | UNIPROT ID: O95466 | |
| Isotype: IgG | Full Name: formin-like 1 | |
| Immunogen Catalog Number: AG13697 | Calculated MW: 122 kDa | |

Applications

Tested Applications:
ELISA

Cited Applications:
WB, IHC

Species Specificity:
human, mouse, rat

Cited Species:
human

Notable Publications

| Author | Pubmed ID | Journal | Application |
|--------|-----------|-----------------|-------------|
| Gui Ma | 37454877 | Mol Cell Probes | IHC, WB |

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data