

For Research Use Only

# NIPA2 Polyclonal antibody

Catalog Number: 20196-1-AP



---

## Basic Information

<b>Catalog Number:</b> 20196-1-AP	<b>GenBank Accession Number:</b> BC011775	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul	<b>GeneID (NCBI):</b> 81614	
<b>Source:</b> Rabbit	<b>Full Name:</b> non imprinted in Prader-Willi/Angelman syndrome 2	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 360 aa, 39 kDa	
<b>Immunogen Catalog Number:</b> AG14128		

---

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human, mouse, rat

---

## Background Information

---

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

---

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data