

For Research Use Only

POTEB-Specific Polyclonal antibody

Catalog Number: 20205-1-AP



Basic Information

Catalog Number:

20205-1-AP

Size:

150ul , Concentration: 253 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_207355

GeneID (NCBI):

339010

UNIPROT ID:

Q6S5H4

Full Name:

POTE ankyrin domain family, member B

Calculated MW:

66 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

POTEB, also named as A26B1 and POTE15, belongs to the POTE family. POTE has an important signaling function in the reproductive system. This antibody is specific to POTEB.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

