

For Research Use Only

# SLIT3-Specific Polyclonal antibody

Catalog Number: 20218-1-AP



## Basic Information

<b>Catalog Number:</b> 20218-1-AP	<b>GenBank Accession Number:</b> NM_003062	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 327 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 6586	
<b>Source:</b> Rabbit	<b>Full Name:</b> slit homolog 3 (Drosophila)	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 168 kDa	

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human

## Background Information

SLIT3, also named as KIAA0814, MEGF5 and SLIL2, may act as molecular guidance cue in cellular migration. It may be mediated by interaction with roundabout homolog receptors. The antibody is specific to SLIT3.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data