

For Research Use Only

# SLIT3-Specific Polyclonal antibody, PBS Only

Catalog Number: 20218-1-PBS



## Basic Information

Catalog Number:

20218-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_003062

GeneID (NCBI):

6586

UNIPROT ID:

O75094

Full Name:

slit homolog 3 (Drosophila)

Calculated MW:

168 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

SLIT3, also named as KIAA0814, MEGF5 and SLIL2, may act as molecular guidance cue in cellular migration. It may be mediated by interaction with roundabout homolog receptors. The antibody is specific to SLIT3.

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20218-1-AP (SLIT3-Specific antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20218-1-PBS in a different storage buffer formulation.