

For Research Use Only

# FBXL5 Polyclonal antibody

Catalog Number: 20295-1-AP



## Basic Information

Catalog Number: 20295-1-AP	GenBank Accession Number: BC030656	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 26234	
Source: Rabbit	Full Name: F-box and Leucine-rich repeat protein	
Isotype: IgG	5	
Immunogen Catalog Number: AG14101	Calculated MW: 691 aa, 79 kDa	

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data