

TSPAN4 Polyclonal antibody

Catalog Number: 20318-1-AP

Basic Information

Catalog Number:

20318-1-AP

Size:

150ul , Concentration: 213 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14172

GenBank Accession Number:

BC000389

GeneID (NCBI):

7106

UNIPROT ID:

O14817

Full Name:

tetraspanin 4

Calculated MW:

238 aa, 26 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

TSPAN4 (also known as NAG2 or TM4SF7) is a membrane protein of the tetraspanin superfamily, which are characterized by the presence of four conserved transmembrane regions. Tetraspanins have been involved in diverse processes such as cell activation and proliferation, adhesion and motility, differentiation, and cancer. TSPAN4 forms complexes with integrins and other tetraspanins (PMID: 9360996).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

