

For Research Use Only

MPDU1 Polyclonal antibody

Catalog Number: 20321-1-AP



Basic Information

Catalog Number: 20321-1-AP	GenBank Accession Number: BC001898	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 340 ug/ml by Nanodrop;	GeneID (NCBI): 9526	Recommended Dilutions: WB: 1:1000-1:4000
Source: Rabbit	UNIPROT ID: O75352	
Isotype: IgG	Full Name: mannose-P-dolichol utilization defect 1	
Immunogen Catalog Number: AG14178	Calculated MW: 247 aa, 27 kDa	
	Observed MW: 20 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human, mouse, rat

Positive Controls:

WB : human placenta tissue, human testis tissue, mouse testis tissue, rat testis tissue

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20321-1-AP (MPDU1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.