

For Research Use Only

MPDU1 Polyclonal antibody, PBS Only

Catalog Number: 20321-1-PBS



Basic Information

Catalog Number:

20321-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14178

GenBank Accession Number:

BC001898

GeneID (NCBI):

9526

UNIPROT ID:

O75352

Full Name:

mannose-P-dolichol utilization defect 1

Calculated MW:

247 aa, 27 kDa

Observed MW:

20 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20321-1-AP (MPDU1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20321-1-PBS in a different storage buffer formulation.