

For Research Use Only

# OSM Polyclonal antibody

Catalog Number: 20345-1-AP

1 Publications



## Basic Information

Catalog Number:

20345-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14152

GenBank Accession Number:

BC011589

GeneID (NCBI):

5008

UNIPROT ID:

P13725

Full Name:

oncostatin M

Calculated MW:

252 aa, 28 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

IHC

Species Specificity:

human

Cited Species:

human

## Notable Publications

Author	Pubmed ID	Journal	Application
Sen Wang	28745433	Cancer Med	IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

