

For Research Use Only

C11orf67 Polyclonal antibody, PBS Only

Catalog Number: 20414-1-PBS

Featured Product



Basic Information

Catalog Number:

20414-1-PBS

Size:

100ug, Concentration: 1mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14258

GenBank Accession Number:

BC002752

GeneID (NCBI):

28971

UNIPROT ID:

Q9H7C9

Full Name:

chromosome 11 open reading frame

67

Calculated MW:

122 aa, 13 kDa

Observed MW:

13 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, ELISA

Species Specificity:

human, mouse, rat

Background Information

C11orf67 is predicted to be involved in positive regulation of fat cell differentiation.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human ovary tumor tissue slide using 20414-1-AP (C11orf67 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 20414-1-PBS in a different storage buffer formulation.



mouse testis tissue were subjected to SDS PAGE followed by western blot with 20414-1-AP (C11orf67 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20414-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human ovary tumor tissue slide using 20414-1-AP (C11orf67 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 20414-1-PBS in a different storage buffer formulation.



WB result of C11orf67 antibody (20414-1-AP; 1:600; incubated at room temperature for 1.5 hours) with sh-Control and sh-C11orf67 transfected HepG2 cells. This data was developed using the same antibody clone with 20414-1-PBS in a different storage buffer formulation.