

For Research Use Only

KBTBD4 Polyclonal antibody

Catalog Number: 20428-1-AP



Basic Information

Catalog Number:

20428-1-AP

Size:

150ul , Concentration: 350 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14261

GenBank Accession Number:

BC002736

GeneID (NCBI):

55709

UNIPROT ID:

Q9NVX7

Full Name:

kelch repeat and BTB (POZ) domain containing 4

Calculated MW:

518 aa, 58 kDa

Observed MW:

50-58 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : A431 cells, C6 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20428-1-AP (KBTBD4 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.