

For Research Use Only

# TMEM101 Polyclonal antibody

Catalog Number: 20473-1-AP



## Basic Information

Catalog Number:  
20473-1-AP

Size:  
150ul

Source:  
Rabbit

Isotype:  
IgG

Immunogen Catalog Number:  
AG14323

GenBank Accession Number:  
BC007438

GeneID (NCBI):  
84336

Full Name:  
transmembrane protein 101

Calculated MW:  
257 aa, 29 kDa

Purification Method:  
Antigen affinity purification

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data