

For Research Use Only

NHE1 Polyclonal antibody

Catalog Number: 20474-1-AP **1 Publications**



Basic Information

| | | |
|--------------------------------------|---|---|
| Catalog Number: 20474-1-AP | GenBank Accession Number: BC012121 | Purification Method: Antigen affinity purification |
| Size: 150ul | GeneID (NCBI): 6548 | |
| Source: Rabbit | Full Name: solute carrier family 9 (sodium/hydrogen exchanger), member 1 | |
| Isotype: IgG | Calculated MW: 815 aa, 91 kDa | |
| Immunogen Catalog Number: AG14324 | | |

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|---------------|-----------|-------------|-------------|
| Shimeng Zhang | 26397063 | Int J Oncol | WB |

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data