

For Research Use Only

NHE1 Polyclonal antibody

Catalog Number: 20474-1-AP

1 Publications



Basic Information

Catalog Number:

20474-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14324

GenBank Accession Number:

BC012121

GeneID (NCBI):

6548

UNIPROT ID:

P19634

Full Name:

solute carrier family 9
(sodium/hydrogen exchanger),
member 1

Calculated MW:

815 aa, 91 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Notable Publications

Author	Pubmed ID	Journal	Application
Shimeng Zhang	26397063	Int J Oncol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

