

For Research Use Only

SNRNP27 Polyclonal antibody, PBS Only

Catalog Number: 20505-1-PBS



Basic Information

Catalog Number:

20505-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14308

GenBank Accession Number:

BC017890

GeneID (NCBI):

11017

UNIPROT ID:

Q8WVK2

Full Name:

small nuclear ribonucleoprotein 27kDa (U4/U6.U5)

Calculated MW:

155 aa, 19 kDa

Observed MW:

23-27 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

SNRNP27 also named U4/U6.U5 small nuclear ribonucleoprotein 27 kDa protein. It is located in the cell nucleus. SNRNP27 may play a role in mRNA splicing.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20505-1-AP (SNRNP27 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20505-1-PBS in a different storage buffer formulation.