

For Research Use Only

ZC4H2 Polyclonal antibody, PBS Only

Catalog Number: 20526-1-PBS



Basic Information

Catalog Number:

20526-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14388

GenBank Accession Number:

BC004411

GeneID (NCBI):

55906

UNIPROT ID:

Q9NQZ6

Full Name:

zinc finger, C4H2 domain containing

Calculated MW:

224 aa, 26 kDa

Observed MW:

26 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20526-1-AP (ZC4H2 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20526-1-PBS in a different storage buffer formulation.