

For Research Use Only

HIST1H3D Polyclonal antibody, PBS Only

Catalog Number: 20532-1-PBS



Basic Information

Catalog Number:

20532-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14394

GenBank Accession Number:

BC031333

GeneID (NCBI):

8351

UNIPROT ID:

P68431

Full Name:

histone cluster 1, H3d

Calculated MW:

136 aa, 15 kDa

Observed MW:

15-25 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 20532-1-AP (HIST1H3D Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20532-1-PBS in a different storage buffer formulation.