

For Research Use Only

FAM111B Polyclonal antibody

Catalog Number: 20552-1-AP

Featured Product



Basic Information

Catalog Number:

20552-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_198947

GeneID (NCBI):

374393

UNIPROT ID:

Q6SJ93

Full Name:

family with sequence similarity 111,
member B

Calculated MW:

85 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

FAM111B, also named as CANP, belongs to the FAM111 family. The antibody is specific to FAM111B. This protein has some isoforms with MW 135 kDa (EAW73833, 1200aa), 85 kDa and 81 kDa. FAM111B is a direct target of p53. It is involved in tendon contracture, myopathy, pulmonary fibrosis and lung adenocarcinoma. (PMID: 26495788, 31114230, 24268661)

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

