

For Research Use Only

FAM111B Polyclonal antibody

Catalog Number: 20552-1-AP

Featured Product



Basic Information

| | | |
|--------------------------------------|---|--|
| Catalog Number: 20552-1-AP | GenBank Accession Number: NM_198947 | Purification Method: Antigen affinity purification |
| Size: 150ul | GeneID (NCBI): 374393 | |
| Source: Rabbit | Full Name: family with sequence similarity 111, member B | |
| Isotype: IgG | Calculated MW: 85 kDa | |

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

FAM111B, also named as CANP, belongs to the FAM111 family. The antibody is specific to FAM111B. This protein has some isoforms with MW 135 kDa (EAW73833, 1200aa), 85 kDa and 81 kDa. FAM111B is a direct target of p53. It is involved in tendon contracture, myopathy, pulmonary fibrosis and lung adenocarcinoma. (PMID: 26495788, 31114230, 24268661)

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data