

For Research Use Only

IKK Alpha Polyclonal antibody

Catalog Number: 20565-1-AP

3 Publications



Basic Information

Catalog Number:

20565-1-AP

Size:

150ul , Concentration: 193 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001278

GeneID (NCBI):

1147

UNIPROT ID:

O15111

Full Name:

conserved helix-loop-helix
ubiquitous kinase

Calculated MW:

85 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human, mouse

Notable Publications

| Author | Pubmed ID | Journal | Application |
|----------------|-----------|---------------------|-------------|
| Rui Wang | 33831690 | Phytomedicine | WB |
| X U Baogui | 36083478 | J Tradit Chin Med | WB |
| Xiao-Long Wang | 34517275 | Int Immunopharmacol | WB |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

