

For Research Use Only

# IKK Alpha Polyclonal antibody

Catalog Number: 20565-1-AP **2 Publications**



## Basic Information

|  |   |  |
|--|---|--|
| <b>Catalog Number:</b><br>20565-1-AP   | <b>GenBank Accession Number:</b><br>NM_001278                     | <b>Purification Method:</b><br>Antigen affinity purification |
| <b>Size:</b><br>150ul , Concentration: 193 µg/ml by Bradford method using BSA as the standard; | <b>GeneID (NCBI):</b><br>1147                                     |  |
| <b>Source:</b><br>Rabbit   | <b>Full Name:</b><br>conserved helix-loop-helix ubiquitous kinase |  |
| <b>Isotype:</b><br>IgG   | <b>Calculated MW:</b><br>85 kDa                                   |  |

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
WB

**Species Specificity:**  
human

**Cited Species:**  
human

## Background Information

### Notable Publications

| Author         | Pubmed ID | Journal             | Application |
|----------------|-----------|---------------------|-------------|
| Rui Wang       | 33831690  | Phytomedicine       | WB          |
| Xiao-Long Wang | 34517275  | Int Immunopharmacol | WB          |

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data