

For Research Use Only

DNMT3B Polyclonal antibody

Catalog Number: 20694-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

20694-1-AP

Size:

150ul , Concentration: 220 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_006892

GeneID (NCBI):

1789

UNIPROT ID:

Q9UBC3

Full Name:

DNA (cytosine-5-)-methyltransferase

3 beta

Calculated MW:

96 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Notable Publications

Author	Pubmed ID	Journal	Application
Sheng Zeng	27720955	Biochim Biophys Acta	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

