

For Research Use Only

DNMT3B Polyclonal antibody

Catalog Number: 20694-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number: 20694-1-AP	GenBank Accession Number: NM_006892	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 220 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 1789	
Source: Rabbit	Full Name: DNA (cytosine-5-)-methyltransferase 3 beta	
Isotype: IgG	Calculated MW: 96 kDa	

Applications

Tested Applications:
WB, ELISA

Cited Applications:
WB

Species Specificity:
human

Cited Species:
human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Sheng Zeng	27720955	Biochim Biophys Acta	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data