

For Research Use Only

EVA1B Polyclonal antibody, PBS Only

Catalog Number:20793-1-PBS



Basic Information

Catalog Number: 20793-1-PBS	GenBank Accession Number: BC006241	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 55194	
Source: Rabbit	UNIPROT ID: Q9NVM1	
Isotype: IgG	Full Name: family with sequence similarity 176, member B	
Immunogen Catalog Number: AG14578	Calculated MW: 165 aa, 18 kDa	
	Observed MW: 25-26 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

EVA1B, also known as FAM176B or C1orf78, is a single-pass membrane protein that belongs to the FAM176 family. EVA1B is an important paralog of the EVA1A and is associated with the tumor immune microenvironment and clinical prognosis

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue was subjected to SDS PAGE followed by western blot with 20793-1-AP (EVA1B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20793-1-PBS in a different storage buffer formulation.