

For Research Use Only

ZKSCAN3 Polyclonal antibody

Catalog Number: 20800-1-AP

Featured Product

3 Publications



Basic Information

Catalog Number:

20800-1-AP

Size:

150ul, Concentration: 347 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14715

GenBank Accession Number:

BC006118

GeneID (NCBI):

80317

UNIPROT ID:

Q9BRRO

Full Name:

zinc finger with KRAB and SCAN domains 3

Calculated MW:

538 aa, 61 kDa

Observed MW:

61-65 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF, IHC

Species Specificity:

human, mouse

Cited Species:

mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Martijn J C van der Lienden	33802460	Int J Mol Sci	IHC
Haihui Pan	28581889	Autophagy	WB
Albert Torra	29628303	Mol Ther	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

