

For Research Use Only

B3GALT4 Polyclonal antibody

Catalog Number: 20817-1-AP



Basic Information

| | | |
|--|---|--|
| Catalog Number: 20817-1-AP | GenBank Accession Number: BC032574 | Purification Method: Antigen affinity purification |
| Size: 150ul , Concentration: 147 µg/ml by Bradford method using BSA as the standard; | GeneID (NCBI): 8705 | |
| Source: Rabbit | Full Name: UDP-Gal:betaGlcNAc beta 1,3-galactosyltransferase, polypeptide 4 | |
| Isotype: IgG | Calculated MW: 378 aa, 42 kDa | |
| Immunogen Catalog Number: AG4603 | | |

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data