

For Research Use Only

MICU2 Polyclonal antibody, PBS Only

Catalog Number:20859-1-PBS



Basic Information

Catalog Number:

20859-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14830

GenBank Accession Number:

BC034965

GeneID (NCBI):

221154

UNIPROT ID:

Q8IYU8

Full Name:

EF-hand domain family, member A1

Calculated MW:

434 aa, 50 kDa

Observed MW:

~40 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20859-1-AP (EFHA1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20859-1-PBS in a different storage buffer formulation.