

For Research Use Only

KLHDC8B Polyclonal antibody, PBS Only

Catalog Number: 20935-1-PBS



Basic Information

Catalog Number: 20935-1-PBS	GenBank Accession Number: BC039083	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 200942	
Source: Rabbit	UNIPROT ID: Q8IXV7	
Isotype: IgG	Full Name: kelch domain containing 8B	
Immunogen Catalog Number: AG15028	Calculated MW: 354 aa, 38 kDa	
	Observed MW: 38 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20935-1-AP (KLHDC8B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20935-1-PBS in a different storage buffer formulation.