

For Research Use Only

FAM55A Polyclonal antibody, PBS Only

Catalog Number:20996-1-PBS



Basic Information

Catalog Number: 20996-1-PBS	GenBank Accession Number: BC029049	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 120400	
Source: Rabbit	UNIPROT ID: Q8N323	
Isotype: IgG	Full Name: family with sequence similarity 55, member A	
Immunogen Catalog Number: AG15021	Calculated MW: 547 aa, 63 kDa	

Applications

Tested Applications:
IHC, Indirect ELISA

Species Specificity:
human

Background Information

FAM55A, also known as NXPE1, belongs to the NXPE family. FAM55A is a poorly characterized gene implicated in cellular processes such as cell cycle regulation and neuronal development.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded colon cancer slide using 20996-1-AP (FAM55A antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 20996-1-PBS in a different storage buffer formulation.