

For Research Use Only

FAM186B Polyclonal antibody

Catalog Number: 21022-1-AP



Basic Information

Catalog Number:

21022-1-AP

Size:

150ul , Concentration: 1200 ug/ml by Nanodrop and 633 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15284

GenBank Accession Number:

BC035621

GeneID (NCBI):

84070

UNIPROT ID:

Q8IYMO

Full Name:

family with sequence similarity 186, member B

Calculated MW:

893 aa, 104 kDa

Observed MW:

120 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, MCF-7 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 21022-1-AP (FAM186B Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.