

For Research Use Only

ABO Polyclonal antibody

Catalog Number: 21047-1-AP



Basic Information

Catalog Number: 21047-1-AP	GenBank Accession Number: BC069595	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 227 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 28	
Source: Rabbit	Full Name: ABO blood group (transferase A, alpha 1-3-N-acetylgalactosaminyltransferase; transferase B, alpha 1-3-galactosyltransferase)	
Isotype: IgG	Calculated MW: 354 aa, 41 kDa	
Immunogen Catalog Number: AG15389		

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data