

For Research Use Only

TBC1D19 Polyclonal antibody, PBS Only

Catalog Number: 21085-1-PBS



Basic Information

Catalog Number: 21085-1-PBS	GenBank Accession Number: BC031642	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 55296	
Source: Rabbit	UNIPROT ID: Q8N5T2	
Isotype: IgG	Full Name: TBC1 domain family, member 19	
Immunogen Catalog Number: AG15323	Calculated MW: 526 aa, 60 kDa	
	Observed MW: 55 kDa	

Applications

Tested Applications:
WB, IHC, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

The function of TBC1D19 remains largely unknown. TBC1D19 may act as a GTPase-activating protein for Rab family protein(s).

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21085-1-AP (TBC1D19 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21085-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 21085-1-AP (TBC1D19 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 21085-1-PBS in a different storage buffer formulation.