

For Research Use Only

FEZ2 Polyclonal antibody, PBS Only

Catalog Number: 21101-1-PBS



Basic Information

Catalog Number:

21101-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15370

GenBank Accession Number:

BC110386

GeneID (NCBI):

9637

UNIPROT ID:

Q9UHY8

Full Name:

fasciculation and elongation protein
zeta 2 (zygin II)

Calculated MW:

353 aa, 40 kDa

Observed MW:

40 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21101-1-AP (FEZ2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21101-1-PBS in a different storage buffer formulation.



WB result of FEZ2 antibody (21101-1-AP; 1:3000; room temperature for 1.5 hours) with wild-type and FEZ2 knockout HAP1 cells. This data was developed using the same antibody clone with 21101-1-PBS in a different storage buffer formulation.



WB result of FEZ2 antibody (21101-1-AP; 1:3000; incubated at room temperature for 1.5 hours) with sh-Control and sh-FEZ2 transfected HepG2 cells. This data was developed using the same antibody clone with 21101-1-PBS in a different storage buffer formulation.