

For Research Use Only

# TTL5 Polyclonal antibody, PBS Only

Catalog Number: 21103-1-PBS



## Basic Information

<b>Catalog Number:</b> 21103-1-PBS	<b>GenBank Accession Number:</b> BC002766	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 100ug, Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 23093	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q6EMB2	
<b>Isotype:</b> IgG	<b>Full Name:</b> tubulin tyrosine ligase-like family, member 5	
<b>Immunogen Catalog Number:</b> AG15373	<b>Calculated MW:</b> 1281 aa, 144 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA  
**Species Specificity:**  
human

## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data