

For Research Use Only

EBF3 Polyclonal antibody, PBS Only

Catalog Number: 21148-1-PBS



Basic Information

Catalog Number:

21148-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15693

GenBank Accession Number:

BC126130

GeneID (NCBI):

253738

UNIPROT ID:

Q9H4W6

Full Name:

early B-cell factor 3

Calculated MW:

596 aa, 65 kDa

Observed MW:

60-70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 21148-1-AP (EBF3 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21148-1-PBS in a different storage buffer formulation.