

For Research Use Only

GALNT5 Polyclonal antibody

Catalog Number: 21225-1-AP



Basic Information

Catalog Number:

21225-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15674

GenBank Accession Number:

BC142676

GeneID (NCBI):

11227

UNIPROT ID:

Q7Z7M9

Full Name:

UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetylglucosaminyltransferase 5 (GalNAc-T5)

Calculated MW:

940 aa, 106 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

